




H.ALLAH DIN SURGICAL

LAPAROSCOPIC INSTRUMENTS



LAPAROSCOPIC INSTRUMENTS, FORCEPS, DISSECTOR

68-1406
68-1409

Forceps, Dissector MARYLAND, Straight,
Double Action, Cross Serrations
HandleLengthWidth
A330mm5mm
A330mm10mm



68-1460

Forceps, Dissector KELLY, Straight, Double
Action
HandleLengthWidth
A330mm5mm

LAPAROSCOPIC INSTRUMENTS, FORCEPS, GRASPING

68-1400

Forceps, Tumor Grasper, Straight, Double
Action, Blunt End Jaw
HandleLengthWidth
C330mm5mm



See LAPAROSCOPIC INSTRUMENTS, HANDLE REFERENCE GUIDE on page LAP-8 for photo reference of the handle type.

LAPAROSCOPIC INSTRUMENTS, FORCEPS, GRASPING



68-1401 Forceps, Fenestrated Grasping, Straight,
Double Action, Slightly Flat Jaw
HandleLengthWidth
C330mm5mm



68-1402 Forceps, Grasper [MAXI], Straight, Double
Action
HandleLengthWidth
C330mm5mm



68-1403 Forceps, Grasping [JOHAN], Straight, Single
Action, Slightly Flat Jaw
HandleLengthWidth
22mmC330mm5mm
68-1404 38mmC330mm5mm



See LAPAROSCOPIC INSTRUMENTS, HANDLE REFERENCE GUIDE on page LAP-8 for photo reference of the handle type.

LAPAROSCOPIC INSTRUMENTS, FORCEPS, GRASPING

68-1405 Forceps, Grasping, [DEBAKEY MARYLAND],
Straight, Double Action, Atraumatic Teeth
HandleLengthWidth
C330mm5mm



68-1407 Forceps, Grasping, Endo Clinch, Straight,
Double Action, Extremely Atraumatic, Serrated
Jaw
HandleLengthWidth
C330mm5mm



68-1408 Forceps, Grasping [ALLIS], Straight, Double
Action
HandleLengthWidth
Rounded JawC330mm5mm
68-1412 Atraumatic,
3X4 TeethC330mm5mm



68-1410 Forceps, Standard Grasping, Straight, Double
Action, Cross Serrated Jaw
HandleLengthWidth
C330mm5mm



68-1411 Forceps, Wave Grasping DG, Straight, Single
Action, Atraumatic, Wave Like Serrated Jaw
HandleLengthWidth
C330mm5mm



See LAPAROSCOPIC INSTRUMENTS, HANDLE REFERENCE GUIDE on page LAP-8 for photo reference of the handle type.

LAPAROSCOPIC INSTRUMENTS, FORCEPS, GRASPING



68-1413
68-1414

Forceps, Grasping [Babcock], Straight, Double
Action
HandleLengthWidth
C330mm10mm
C330mm5mm



68-1415

Forceps, Grasping, TC [Collin], Straight,
Double Action, 30MM Jaw
HandleLengthWidth
C330mm10mm



68-1420

Forceps, Clamp Grasping [Mixer], 90 Degree
Angled, Double Action
HandleLengthWidth
C330mm5mm




68-1421

Forceps, Grasping [Kocher], Straight, Double
Action, Axial Serrations, 60mm Jaw
HandleLengthWidth
C330mm10mm

See LAPAROSCOPIC INSTRUMENTS, HANDLE REFERENCE GUIDE on page LAP-8 for photo reference of the handle type.

LAPAROSCOPIC INSTRUMENTS, FORCEPS, GRASPING, DETACHABLE



68-1416 Forceps, Grasping, Straight
Single Action, Long Atraumatic, Serrated
HandleLengthWidth
B450mm5mm

68-1417 Double Action, Blunt Tip-w/Reservoir, Serrated
HandleLengthWidth
B450mm5mm



68-1418 Forceps, Grasping, [Maxi], Straight, Double
Action
HandleLengthWidth
B450mm5mm



68-1419 Forceps, Grasping, Alligator, Straight, Double
Action
HandleLengthWidth
B450mm5mm



LAPAROSCOPIC INSTRUMENTS, FORCEPS, SPOON

	Forceps, Straight, Single Action	
	HandleLength	Width
68-1461	Biopsy Spoon C330mm	5mm
68-1462	Large SpoonC330mm	10mm



See LAPAROSCOPIC INSTRUMENTS, HANDLE REFERENCE GUIDE on page LAP-8 for photo reference of the handle type.

LAPAROSCOPIC INSTRUMENTS, NEEDLE HOLDER



68-1555



68-1556



68-1557




68-1558



68-1563



68-1564



68-1565

Needleholder, Straight, Single Action		Handle	Length	Width
68-1555	Axial Handles	F	330mm	5mm
68-1556	Axial Handles	F	450mm	5mm
68-1557	DG Handles	D	330mm	5mm
68-1558	DG Handles	D	450mm	5mm
68-1563	Tungsten Carbide Inserts	E	330mm	5mm
68-1564	Left Curved, Tungsten Carbide Inserts	G	330mm	5mm
68-1565	Right Curved, Tungsten Carbide Inserts	G	330mm	5mm

See LAPAROSCOPIC INSTRUMENTS, HANDLE REFERENCE GUIDE on page LAP-8 for photo reference of the handle type.

LAPAROSCOPIC INSTRUMENTS, SCISSORS

Scissors [METZENBAUM], Double Action
Handle LengthWidth
68-1580 Straight, SerratedA330mm5mm
68-1581 Curved, SerratedC330mm5mm
68-1582 Curved, Medium
BladesA330mm5mm



68-1580




68-1582

68-1583 Scissors [MINI-METZENBAUM], Curved
Double Action
HandleLengthWidth
A330mm5mm



LAPAROSCOPIC INSTRUMENTS, TRAY

68-1599 Laparoscopic Instrument Tray. Includes
Laparoscopic Instrument Rack



LAPAROSCOPIC INSTRUMENTS, ACCESSORIES, MAT

68-1598 Mat, Pin For Laparoscopic Instrument Tray
[68-1599]

See LAPAROSCOPIC INSTRUMENTS, HANDLE REFERENCE GUIDE on page LAP-8 for photo reference of the handle type.

LAPAROSCOPIC INSTRUMENTS, HANDLE REFERENCE GUIDE

A Handle, Peek, Rotating, Cleaning Port



E Handle, Stainless Steel, InLine, Ring, Ratchet



B Handle, Peek, Rotating, Ratcheted, Cleaning Port



F Handle, Stainless Steel, Axial, Ratchet



C Handle, Peek, Rotating, Ratcheted, Cleaning Port



G Handle, Stainless Steel, 30° Axial, Ratchet



D Handle, Stainless Steel, Rotating, Slide Ratchet

