



# H.ALLAH DIN SURGICAL

---Laryngoscope Sets---



## Conventional Macintosh set

Code	Fig	Working Length (mm)
6004 A	0	55
6004 B	1	70
6004 C	2	78
6004 D	3	108
6004 E	4	133

Supply in hard case  
Require 2X C Batteries



## Fiber Optic Macintosh set

Code	Fig	Working Length (mm)
6006 A	0	55
6006 B	1	70
6006 C	2	78
6006 D	3	108
6006 E	4	133

Supply in hard case  
Require 2X C Batteries



# H. ALLAH DIN SURGICAL



**Fiber Optic Miller set**

Code	Fig	Working Length (mm)
6007 A	00	45
6007 B	0	55
6007 C	1	80
6007 D	2	130
6007 E	3	170
6007 F	4	180

Supply in hard case  
Require 2X C Batteries



**Fiber Optic Miller set**

Code	Fig	Working Length (mm)
6005 A	00	45
6005 B	0	55
6005 C	1	80
6005 D	2	130
6005 E	3	170
6005 F	4	180

Supply in hard case  
Require 2X C Batteries



# H.ALLAH DIN SURGICAL



**GM Dermatoscope** (German Technology)

**Supply:** with soft pouch

**Require:** 2X C size batteries

## Otoscope Set



### Fiber Optical Otoscope

Bayonet Connected System Fiber optic light halogen, illumination 2.5v  
High quality big swiveling viewing window with 3x magnification lens,

**Supply:** 10 single use mini specula

**Size:** 2.5 & 4mm Supply with soft pouch

**Require:** 2x C Batteries



# H.ALLAH DIN SURGICAL



## Conventional Otoscope

Bayonet Connected System realizable and indestructible high quality Otoscope  
Invaluable for everyday ear throat and nose use regulated easily neutral 2.5v vacuum light  
(German Bulbs) high quality magnification lens & 3 specula  
(**Size:** 2.5, 3.5 & 4.5mm Supply with soft pouch)  
**Require:** 2x C Batteries



## Mini Fiber Optic Otoscope

High quality professional plastic Otoscope with extra bright 2.5v, xenon  
illumination light provided through fiber optic system, high quality swiveling viewing window  
3x magnification lens. Output for pneumatic membrane test.  
**Supply:** 10 single use mini specula  
**Size:** 2.5 & 4mm Supply with soft pouch  
**Require:** 2x AA Batteries



# H.ALLAH DIN SURGICAL



## **Conventional Otoscope**

Bayonet Connected System Realizable and indestructible high quality Otoscope  
Invaluable for everyday ear throat and nose use regulated easily neutral 2.5v vacuum light  
(German Bulbs) high quality magnification lens & 3 specula  
(**Size:** 2.5, 3.5 & 4.5mm Supply with soft pouch)

**Require:** 2x C Batteries



## H.ALLAH DIN SURGICAL



### **Ophthalmoscope Otoscope kit**

Ophthalmoscope 2.5v light for tissue color & consistent,  
long testing illumination 24 color coated lenses 25 to + 40 for excellent resolution  
(Red for - Green for +) Rubber eye brow rest to prevent scratching of eye-glasses  
Illumination lens dial for clear identification of diopter setting.

**Supply:** with parker otoscope

**Packing:** Pouch or Plastic Box

**Require:** 2x C Batteries



## H.ALLAH DIN SURGICAL



### **Mini Ophthalmoscope + Fiber Optic Otoscope Set**

5 Focus Ophthalmoscope (small, large, blue, half moon & net focus)  
2.5v Halogen light for true tissue color & consistence, long testing illumination 16 color coated lenses - 25 + 25 for excellent resolution (Red for - Green for +)  
Rubber eye brow rest to prevent scratching of eye-glasses illumination lens dial for clear identification for adopter setting.

**Supply:** Hard plastic box packing

**Require:** 2x AA batteries



# H.ALLAH DIN SURGICAL



## Small Diagnostic Set

Bayonet Contacting Fitting Regular plastic containing one 2.5v German Bulb Otoscope 3 ear Specula (2.5, 3.5 & 4.5mm) One ophthalmoscope 1 bend arm throat lamp, 2 laryngeals mirrors (Size 4 & 5) Tongue spatula support for single use tongue depressor, adjustable nasal speculum , Battery handle chrome plated rheostat for adjustment of illumination.

**Require:** 2x C Batteries

**Supply:** Storage black plastic case





# H.ALLAH DIN SURGICAL



**Diagnostic Set**

Light weight robust zipper case containing  
1 Otoscope &  
1 Ophthalmoscope  
1 Bend arm throat lamp  
2 Laryngeals mirrors (size 3 & 4)  
1 Tongue spatula  
1 Adjustable Nasal Speculum  
1 Plastic Converter for mini specula fitting  
2 Bottles spear mini specula (2.5mm & 4mm)  
**Require:** 2x C Batteries